

## Photo-physical Properties of Ancient and Modern Artwork Pigments

Gianluca Accorsi<sup>1</sup>, Giovanni Verri<sup>2</sup>, Georgios Arnaoutakis<sup>3</sup>

<sup>1</sup> Italian National Research Council, CNR NANO, Lecce, Italy,

<sup>2</sup> The Courtauld Institute of Art, London, UK, <sup>3</sup> Edinburgh Instruments Ltd, Livingston, UK

AN\_P23 v.2

DATE 7 Dec. 15



**Edinburgh Instruments Ltd**  
2 Bain Square,  
Kirkton Campus,  
Livingston, United Kingdom  
EH54 7DQ

**Telephone (UK)**  
+44(0) 1506 425 300  
**Telephone (North America)**  
1-800-323-6115

**Email**  
[enquiries@edinst.com](mailto:enquiries@edinst.com)  
**Website**  
[www.edinst.com](http://www.edinst.com)

AN\_P23

Copyright © 2015. Edinburgh Instruments Ltd. All rights reserved

1

## Introduction

The interest in non-invasive investigations for the study and identification of painting materials has greatly increased during the past two decades.<sup>1</sup> The nature of museum objects is such that sampling is always kept to a minimum, thus, a spectroscopic approach is ideal for this purpose. In this respect, an ancient and a modern pigment, cuprorivaite,  $\text{CaCuSi}_4\text{O}_{10}$  (Egyptian blue)<sup>2</sup> and barium manganate (VI) sulphate compound (manganese blue)<sup>3</sup> respectively, have been photo-physically investigated. The investigation includes qualitative information, *i.e.* emission and excitation spectra, as well as quantitative that is the photoluminescence quantum yields (PLQY) that have been obtained in the near-infrared (NIR) region.

## Methods & Materials

Excitation and emission spectra were measured using an FLS980 Fluorescence Spectrometer equipped with a 450 W Xe lamp with double excitation and emission monochromators. Near-infrared detectors (Hamamatsu) were used for the detection of both samples. PLQY have been calculated by corrected emission spectra obtained by using a barium sulphate coated integrating sphere, following the procedure described by De Mello *et al.*<sup>4</sup> Experimental uncertainties were estimated to be  $\pm 20\%$  for emission quantum yields,  $\pm 2\text{nm}$  and  $\pm 5\text{nm}$  for absorption and emission peaks, respectively.

## Results - Discussion

Figure 1 reports the excitation and luminescence spectra of the two samples.




Figure 1: Excitation and emission spectra of ancient, Egyptian blue, and modern, manganese blue, pigments. The energy levels corresponding to each transition are also displayed.

Egyptian blue shows two different electronic transitions ( ${}^2B_{1g}-{}^2E_g$  and  ${}^2B_{1g}-{}^2A_{1g}$ ) that can be assigned to  $Cu^{2+}$  ions, which are expected to be the only luminescent components of cuprorivaite.<sup>1</sup> On the other hand, those of Manganese Blue can be attributed to ligand-field (LF) transitions ( ${}^2E-{}^2T_2$ , c.a 800 nm - 900 nm) and to ligand-to-metal charge-transfer bands (between 600 nm and 800 nm) of the  $MnO_4^{2-}$  unit.<sup>5</sup>

The luminescence profiles have peaks at 920 nm and 1300 nm for Egyptian blue and manganese blue, respectively. Consequently, the corresponding PLQYs have been obtained. Manganese blue showed a quantum yield of  $\Phi=0.5\%$ <sup>6</sup> attributed to luminescence quenching factors that occur at lower energy. By contrast, Egyptian blue revealed to be a very strong NIR emitter ( $\Phi=10.5\%$ )<sup>1</sup> that, to the best of our knowledge, has the highest quantum efficiency for a molecule-level chromophore in the 800 nm – 1100 nm range.

## Conclusions

The photo-physical properties of an ancient and a modern art pigment were investigated by means of fluorescence spectroscopy. In addition, the PLQY of the pigments has been measured, as a quantitative means of their photoluminescence properties and demonstrate the effectiveness of this method as a minimally invasive means of measuring samples of cultural significance.

## References

1. Accorsi, G.; Verri, G.; Bolognesi, M.; Armaroli, N.; Clementi, C.; Miliani, C.; Romani, A. The exceptional near-infrared luminescence properties of cuprorivaite (Egyptian blue). *Chem. Commun.* 2009, (**23**), 3392-3394.
2. (a) *Artist's Pigments – A Handbook of Their History and Characteristics*. Elisabeth West Fitzhugh editor: 1997; Vol. 3; (b) Berke, H. The invention of blue and purple pigments in ancient times. *Chem. Soc. Rev.* 2007, **36 (1)**, 15-30.
3. Eastaugh, N.; Walsh, V.; Chaplin, T.; Siddall, R. *Barium Manganate(IV) Sulfate*. Butterworth-Heinemann: Oxford (UK), 2008.
4. De Mello, J. C.; Wittmann, H. F.; Friend, R. H. An improved experimental determination of external photoluminescence quantum efficiency. *Adv. Mater.* 1997, **9 (3)**, 230-236.
5. Chen, Z. L.; Schmalke, H. W.; Fox, T.; Blacque, O.; Berke, H. NH-functionalized tungsten complexes of 2-(dimethylphosphino)imidazole. *J. Organomet. Chem.* 2007, **692 (22)**, 4875-4885.
6. Accorsi, G.; Verri, G.; Acocella, A.; Zerbetto, F.; Lerario, G.; Gigli, G.; Saunders, D.; Billinge, R. Imaging, photophysical properties and DFT calculations of manganese blue (barium manganate(VI) sulphate) - a modern pigment. *Chem. Commun.* 2014, **50 (97)**, 15297-15300.

## **Proprietary notice**

Words and logos marked with ® or ™ are registered trademarks or trademarks owned by EDINBURGH INSTRUMENTS Limited. Other brands and names mentioned herein may be the trademarks of their respective owners.

Neither the whole nor any part of the information contained in, or the product described in, this document may be adapted or reproduced in any material form except with the prior written permission of the copyright holder.

The product described in this document is subject to continuous developments and improvements. All particulars of the product and its use contained in this document are given by EDINBURGH INSTRUMENTS in good faith. However, all warranties implied or expressed, including but not limited to implied warranties of merchantability, or fitness for purpose, are excluded.

This document is intended only to assist the reader in the use of the product. EDINBURGH INSTRUMENTS Limited shall not be liable for any loss or damage arising from the use of any information in this document, or any error or omission in such information, or any incorrect use of the product.

## **Confidentiality status**

This document is Open Access. This document has no restriction on distribution.

## **Feedback on this Application Note**

If you have any comments on this Application Note, please send email to [sales@edinst.com](mailto:sales@edinst.com) giving:

the document title

the document number

the page number(s) to which your comments refer

an explanation of your comments.

General suggestions for additions and improvements are also welcome.