

## Introduction

The study of materials emitting in the mid-infrared (MIR) region, particularly around 3  $\mu\text{m}$ , is of significant interest due to their potential applications as gain media in lasers. Laser's utilising Er<sup>3+</sup>-doped media have diverse applications in soft-tissue medicine as their emission wavelength closely aligns with the fundamental stretching vibration of O-H at 2.95  $\mu\text{m}$ . Previous studies have demonstrated the ability of these lasers to selectively excite or ablate protein-containing tissues, bones, and lipid-rich tissues within this spectral region while minimising collateral damage.<sup>1</sup>

The scope of their applications also extends beyond the medical field. The strong absorption of CO<sub>2</sub> at 2.8  $\mu\text{m}$ , CO at 2.4  $\mu\text{m}$ , and NO<sub>2</sub> at 2.9  $\mu\text{m}$  can be utilised in laser-based systems for effective environmental monitoring.<sup>1</sup> Er-doped glass is also widely used as a fibre amplifier for all fibre communication. Using the photoluminescence peak of Er<sup>3+</sup> at 1.5  $\mu\text{m}$ , Er-doped fibres have been developed as the primary amplification medium for C-band (1.530 – 1.565  $\mu\text{m}$ ) and L-band (1.565 – 1.625  $\mu\text{m}$ ) optical amplifiers.<sup>2</sup>

In this Application Note, the photoluminescence spectrum and photoluminescence lifetime of an Er<sup>3+</sup>-doped ZnF<sub>2</sub>-BaF<sub>2</sub>-SrF<sub>2</sub>YF<sub>3</sub> (ZBSY-e) glass are characterised using an Edinburgh Instruments FLS1000 Photoluminescence Spectrometer.

## Materials and Methods

The sample was a 2.5 cm × 1.5 cm piece of Er<sup>3+</sup> doped fluoride glass. The photoluminescence properties were characterised using an Edinburgh Instruments FLS1000 equipped with a 2 W 980 nm laser diode with a pulse modulation (PM-2) box as an excitation source. The detector was a dual-mode MIR InAs-3100 which has a spectral range of 1.2 – 3.1  $\mu\text{m}$ . The glass sample was held using the N-J01 Front Face Sample Holder of the FLS1000.



Figure 1 Edinburgh Instruments FLS1000 Photoluminescence Spectrometer.

## NIR & MIR Spectra

First, the photoluminescence spectrum of the glass was acquired using the 980 nm laser source in CW mode for excitation. Intense emission peaks were observed at around 1.55  $\mu\text{m}$  and 2.75  $\mu\text{m}$  (Figure 2).




Figure 2 Emission spectrum of the Er<sup>3+</sup> doped ZBSY-e fluoride glass. The noise of 2.75  $\mu\text{m}$  peak is from atmospheric MIR absorption lines.

These peaks are assigned to the  $4I_{13/2} \rightarrow 4I_{15/2}$  and the  $4I_{11/2} \rightarrow 4I_{13/2}$  transitions, respectively, as shown in Figure 3. In this context, CR denotes the non-radiative process of cross-relaxation.




Figure 3 Energy-transfer process between excited Er<sup>3+</sup> ions.

The spectrum agrees with previous reports, which examined the efficiency of laser transitions of the rare-earth cations. It is suggested that with optimal engineering of the Er<sup>3+</sup> energy levels and the fibre laser resonator, a high Er<sup>3+</sup> concentration and effective cooling of the host fluoride glass, the commercialisation of high-power MIR lasers can be achieved.<sup>1</sup>

## NIR & MIR Time-Resolved Decays

Next, the photoluminescence decays were measured using the 980 nm laser in pulsed mode at a repetition rate of 10 Hz. The laser pulse width was adjusted using the PM-2 box and was set to a pulse width of 1 ms. The photoluminescence decays of the  $^4I_{1/2} \rightarrow ^4I_{13/2}$  transition at 2.75  $\mu\text{m}$  and  $^4I_{13/2} \rightarrow ^4I_{15/2}$  transition at 1.55  $\mu\text{m}$  are shown in Figure 4. The 2.75  $\mu\text{m}$  decay was fitted using a single exponential model, yielding a lifetime of 7.08 ms. The 1.55  $\mu\text{m}$  decay has a rising component in addition to the decay and was therefore fitted using a double exponential model (one rising component and one decay component) which gave a decay lifetime of 14.2 ms.




Figure 4 Photoluminescence decays of the Er<sup>3+</sup>-doped ZBSY-e fluoride glass at 1.55  $\mu\text{m}$  and 2.75  $\mu\text{m}$ .

## Conclusion

This application note demonstrates the capability of the FLS1000 for NIR and MIR photoluminescence spectra and lifetime measurements. The spectra and lifetimes of an Er<sup>3+</sup> doped fluoride glass were measured, which is an important characterisation step when developing rare-earth-based applications, including lasers, medicine, environmental monitoring, and telecommunications.

## References

1. S. D. Jackson. Towards high-power mid-infrared emission from a fibre laser. *Nat. Photonics* 6, 423–431 (2012).
2. Q. Wang, et al. Enhancement of lifetime in Er-doped silica optical fiber by doping Yb ions via atomic layer deposition. *Opt. Mater. Express* 10, 397 (2020).



For more information, please contact:

+44 (0) 1506 425 300  
[sales@edinst.com](mailto:sales@edinst.com)  
[www.edinst.com](http://www.edinst.com)